



**Presto Stantest Pvt. Ltd.**

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India  
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

• Faridabad • Sonapat • Kolkata • Mumbai • Pune • Ahmedabad  
• Chennai • Bangalore • Hyderabad

[www.prestogroup.com](http://www.prestogroup.com)



## Carbon Black Dispersion Test Apparatus



### Model CBD - 501

**Presto Carbon Black dispersion test apparatus** is mainly used for the determination of carbon black dispersion. It finds wide application in the plastic industry. The instrument comprises of a hot plate which is 8 inch in diameter. The temperature range can be adjusted between ambient and 200°C. Presto carbon black dispersion test apparatus is fitted with a 6 inch size microscope with a magnification power of 100 x and 450 x. The instrument is manufactured as per IS 4984 and IS 4985.





# Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India  
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Sonapat • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad

[www.prestogroup.com](http://www.prestogroup.com)



## Key Specifications:

|                        |                                |
|------------------------|--------------------------------|
| Size of hot plate      | 8 inch diameter.               |
| Temperature range      | Upto 200°C                     |
| Temperature Controller | Thermostate                    |
| Magnification          | 100 X, 450 X. (Any one select) |
| Power                  | 230V, 50 Hz, Single phase.     |

## Highlights:

|                            |  |
|----------------------------|--|
| Finish                     | Powder coated Havel Gray & Blue combination finish and bright chrome / zinc plating for corrosion resistant finish |
| Size of Micro Scope Screen | 6 Inch Diameter  |

Thank you customers for choosing us as your partners in growth !



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

## OUR OTHER PRODUCTS :-

|                        |                         |                        |                    |                    |                     |
|------------------------|-------------------------|------------------------|--------------------|--------------------|---------------------|
|                        |                         |                        |                    |                    |                     |
| Box Compression Tester | Scuff Resistance Tester | Tensile Tester Digital | Spectrocolorimeter | Cobb Sizing Tester | Digital GSM Balance |

