





## Cobb Sizing Tester




### Model No. PCS - 23

**Cobb Sizing Tester** is a consistent tool useful for checking the water absorptiveness (Cobb value) of a given test sample. Used to check the quantity of water absorbed by various sample like paper, cardboard, corrugated sample in a given time rapidly and accurately. PRESTO Cobb Sizing Tester is designed as per standards IS 1060, TAPPI –T441.


- High Quality rugged structure with corrosion resistant main body
- Light weight and resilient
- Fast Locking Mechanism
- Single handed hassle free operation
- Highly accurate test results
- High Grade stainless steel test area



Top View



Front View



Side View



# Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India  
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Sonipat • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad

[www.prestogroup.com](http://www.prestogroup.com)



## Key Specifications:

Sample Diameter	130 mm (approx)
Depth of cylinder	60 mm
Diameter of rubber pad	145 mm
Water level indicator in cylinder	10 mm and 25 mm
Weight of roller	10Kg (approx)

## Highlights:

No. of Roller	1
Locking screw	Yes
Material of Main Body	Stainless Steel / mild steel

Used in various Industries like paper and packaging, textile, films, foils etc.

## Optional Features:

Combine with Digital Weighing balance for a complete Kit.

Thank you customers for choosing us as your partners in growth !



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

## OUR OTHER PRODUCTS :-

